

WHEELCHAIR
Platform

WITH CORRECTA UNIT



DENTAL EQUIPMENTS

*Confidential, Real quality
Since 1943*

WHEELCHAIR Platform

WITH CORRECTA UNIT

On dentist side, control panel, which is produced by hygienic material, provides activating cuspidor, cup filler, operating light and commanding also platform positioning.


TEKMİL SANAYİ MÜESSESESİ

Certified by :



DIN EN ISO 9001:2008
DIN EN ISO 13485:2003
Medical Devices Directives 93-42

CE 1984



With security switches prevents foot from getting stuck.



Multi Axis Planet headrest with free movement capability allows the most convenient position for patients.

Esthetic and ergonomic Correcta dentist tablet has 5 exits. Each exit can be designed as dentist requested instrument (scaler, air/elec. Micromotor, fiber optic etc). Each exit's air and water regulations are separate.

Dental Lights optional photocell or, light emitting diode (LED) and movable head systems



OPTIONAL

Peristaltic Pump for surgical operations with serum holder.

OPTIONAL

Speed and direction information on the LCD display



Aspirator's filter can be removed easily for cleaning



With one button, it is possible to have water directly from unit distilled water tank or from city water



Removable and easy cleanable magnetic instruments pillow.



Silicon cover can be removed and completely autoclavable.



DENTAL EQUIPMENTS

TEKMİL TIBBİ ARAÇ ve GEREÇLERİ TİCARET ve SANAYİ LTD. ŞTİ.

İnönü Mahallesi Çayır Sokak No. 43 - Harbiye - İstanbul / TURKEY

Phone : +90 212 233 78 37

Telefax : +90 212 231 72 77

e-mail : info@tekmil.com.tr

230 49 86

233 42 09

web : www.tekmil.com.tr