



Confidential, real quality since 1943

Epica

SURGERY BOX



**DENTAL
EQUIPMENT**




From Diverso Box control panel;
possible to increase/decrease and
save 3 program of the micromotor speed
and rotate right /left.



Silicone cover can be easily
removed and cleaned,
completely autoclavable

Esthetic and ergonomic
Mobil Box has 5 exits.
Each exit designed as the
instrument desired.

Each exit has air and water
regulator commanded separately.
Chair movements control is optional.



The hoses could be detachable
and autoclavable for more
hygienic working conditions



Water tank removed easily
from the unit part
for filling and cleaning

For surgical operation
it is possible to connect
peristaltic pump with
implant micromotor
& implant contra angle



Pneumatic pedal
with the chip blower system
and water spray command)



Large Tray with the silicone cover



Support arm with 2 exits aspirator surgical, possibility to prepare the other 2 exits,
according to dentist demand.
Support arm produced from a material can be cleaned easily.

Silicone cover can be easily
removed and cleaned,
completely autoclavable



The cuspidor of Epica Surgery Box
is completely ceramic,
hygienic and cleaned easily,
can move right and left side



Rotation of the support arm with 4 exits,
provides a working comfortable to the dentist.
The movements of the chair, cuspidor, and cup
controlled from the control panel.



The exits of cuspidor and glass
can be removed for cleaning



Aspirator's hose and filter
can removed easily for cleaning



Serum connection arm can be extended as needed.

Faro Alya



Faro Maia



Faro Edi



G.Comm Vision




Halogen, Led, photocell and 3rd axis movable head system options




Movable body of unit

Epica
SURGERY BOX





Ergonomic and esthetic Elegant Chair,
with the backrest and the seat movable,
provides to the dentist a working position as desired.



To the chair's movements
regulated by program or
manually, there is also
trendelenburg position
with hip movement

Electric - Pneumatic pedal
with command
of the movements of chair
(with the chip blower system)



Pneumatic pedals
(with the chip blower system and
water spray command)



Soft upholstery option



Narrow and wide
backrest option



Command fixed to the unit's floor,
provides to the dentist,
the possibility to change
the movements of the chair
with foot control, during working
in the mouth of the patient.



Possibility to add movable right arm



Doctor stool with foot ring
and chromium base




Dental stool
(saddle-shaped seat)



The backrest is movable:
up- down, forward - back.
Fixed to the position
as desired



Multifunction
doctor stool
with foot ring and
chromium base



DIN EN ISO 9001
DIN EN ISO 13485
Medical Devices Directives 93-42

CE 1984



DIŞ HEKİMLİĞİ CİHAZLARI

İnönü Mahallesi Çayır Sokak No.43 Harbiye-Şişli-İstanbul/ TURKEY

+90 212 233 7837

230 4986

+90 212 231 7277

+90 545 96367 26

satis@tekmil.com.tr

www.tekmil.com.tr